

Product Datasheet

IFT57 antibody (orb627939)

Description	Rabbit polyclonal antibody to IFT57
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	intraflagellar transport 57 homolog(Chlamydomonas)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	49kd
Uniprot ID	Q9NWB7
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20- 1:200
Expiration Date	12 months from date of receipt.