

Product Datasheet

IFT88 antibody (orb627943)

Description	Rabbit polyclonal antibody to IFT88
Species/Host	Rabbit
Reactivity	Canine, Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	intraflagellar transport 88 homolog(Chlamydomonas)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	88-95 kDa
Uniprot ID	Q13099
Dilution Range	WB: 1:200-1:2000, IP: 1:200-1:2000, IHC-P: 1:20-1:200, IF: 1:20-1:200
Expiration Date	12 months from date of receipt.