

Product Datasheet

IGF2BP2 antibody (orb627946)

Description	Rabbit polyclonal antibody to IGF2BP2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	insulin-like growth factor 2 mRNA binding protein 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	66 kDa
Uniprot ID	Q9Y6M1
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:50-1:100
Expiration Date	12 months from date of receipt.