

Product Datasheet

IL1 Beta antibody (orb627980)

Description	Mouse monoclonal antibody to IL1 Beta
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	interleukin 1, beta
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	35kDa
Uniprot ID	P01584
Dilution Range	WB: 1:1000-1:4000, IHC-P: 1:50-1:500
Expiration Date	12 months from date of receipt.