

Product Datasheet

IL-17 antibody (orb627995)

Description	Mouse monoclonal antibody to IL-17
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, FC, WB
Immunogen	interleukin 17A
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	32-35 kDa
Uniprot ID	Q16552
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.