

Product Datasheet

IL-27 antibody (orb628020)

Description	Rabbit polyclonal antibody to IL-27
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	interleukin 27
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	31 kDa
Uniprot ID	Q8NEV9
Dilution Range	WB: 1:2000-1:10000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.