

Product Datasheet

INO80 antibody (orb628066)

Description	Rabbit polyclonal antibody to INO80
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	INO80 homolog(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	177,150 kDa
Uniprot ID	Q9ULG1
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:2000
Expiration Date	12 months from date of receipt.