

Product Datasheet

KCNF1 antibody (orb628202)

Description	Rabbit polyclonal antibody to KCNF1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	potassium voltage-gated channel, subfamily F, member 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q9H3M0
Dilution Range	IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.