

Product Datasheet

KIF26B antibody (orb628259)

Description	Rabbit polyclonal antibody to KIF26B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	kinesin family member 26B
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	224-230 kDa or 17-22 kDa
Uniprot ID	Q2KJY2
Dilution Range	WB: 1:500 - 1:1000; IHC: 1:50 - 1:200
Expiration Date	12 months from date of receipt.