

Product Datasheet

LGALS3BP antibody (orb628424)

Description	Rabbit polyclonal antibody to LGALS3BP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	lectin, galactoside-binding, soluble, 3 binding protein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	65-70 kDa
Uniprot ID	Q08380
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:2000
Expiration Date	12 months from date of receipt.