

Product Datasheet

MED10 antibody (orb628720)

Description	Rabbit polyclonal antibody to MED10
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	mediator complex subunit 10
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	16 kDa
Uniprot ID	Q9BTT4
Dilution Range	WB: 1:200 - 1:1000
Expiration Date	12 months from date of receipt.