



Product Datasheet

Stathmin (phospho-Ser16) antibody (orb6288)

Description	Stathmin (phospho-Ser16) antibody
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	IHC-P
Immunogen	KLH conjugated synthesised phosphopeptide derived from human Stathmin 1 around the phosphorylation site of Ser16 (RA(p-S)GQ)
Target	Phospho-STMN1 (Ser16)
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20°C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	16 kDa
Uniprot ID	P16949
Dilution Range	IHC-P=1:100-500, WB=1:500-2000, IHC-F=1:100-500, IF=1:100-500, ELISA=1:5000-10000
Expiration Date	12 months from date of receipt.



IHC-P analysis of human brain glioma tissue using Stathmin (phospho-Ser16) antibody (orb6288).



IHC-P analysis of mouse ovary tissue using Stathmin (phospho-Ser16) antibody (orb6288).