

Product Datasheet

MKI67IP antibody (orb628822)

Description	Rabbit polyclonal antibody to MKI67IP
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	MKI67 FHA domain interacting nucleolar phosphoprotein
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	34kd
Uniprot ID	Q9BYG3
Dilution Range	WB: 1:1000-1:4000
Expiration Date	12 months from date of receipt.