

Product Datasheet

MRPL28 antibody (orb628942)

Description	Rabbit polyclonal antibody to MRPL28
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	mitochondrial ribosomal protein L28
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30 kDa
Uniprot ID	Q13084
Dilution Range	WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50- 1:100, IP: 1:50-1:200
Expiration Date	12 months from date of receipt.