

Product Datasheet

MUC16 antibody (orb629043)

Description	Mouse monoclonal antibody to MUC16
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	mucin 16, cell surface associated
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	90-100 kDa
Uniprot ID	Q8WXI7
Dilution Range	IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.