

Product Datasheet

MYC tag antibody (orb629066)

Description	Mouse monoclonal antibody to MYC tag
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	MYC-tag
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
Dilution Range	WB: 1:2000-1:20000, IP: 1:2000-1:20000, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.