

Product Datasheet

MYC-tag antibody (orb629071)

Description	Rabbit polyclonal antibody to MYC-tag
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IP, WB
Immunogen	MYC-tag
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Dilution Range	WB: 1:1000-1:8000, IP: 1:500-1:2000 , IF: 1:10-1:100
Expiration Date	12 months from date of receipt.