

Product Datasheet

MYH9 antibody (orb629085)

Description	Mouse monoclonal antibody to MYH9
Species/Host	Mouse
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	myosin, heavy chain 9, non-muscle
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	224 kDa
Uniprot ID	P35579
Dilution Range	WB: 1:1000-1:5000
Expiration Date	12 months from date of receipt.