

Product Datasheet

NEIL1 antibody (orb629234)

Description	Rabbit polyclonal antibody to NEIL1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	nei endonuclease VIII-like 1(E. coli)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	44kd
Uniprot ID	Q96FI4
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.