

Product Datasheet

NETO2 antibody (orb629247)

Description	Rabbit polyclonal antibody to NETO2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	neuropilin(NRP) and tolloid(TLL)-like 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	59kd
Uniprot ID	Q8NC67
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.