

## Product Datasheet

### NEURL antibody (orb629250)

|                            |                                                                                                                     |
|----------------------------|---------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>         | Rabbit polyclonal antibody to NEURL                                                                                 |
| <b>Species/Host</b>        | Rabbit                                                                                                              |
| <b>Reactivity</b>          | Human, Mouse, Rat                                                                                                   |
| <b>Conjugation</b>         | Unconjugated                                                                                                        |
| <b>Tested Applications</b> | ELISA, IP, WB                                                                                                       |
| <b>Immunogen</b>           | neuralized homolog(Drosophila)                                                                                      |
| <b>Form/Appearance</b>     | liquid                                                                                                              |
| <b>Storage</b>             | PBS with 0.02% sodium azide and 50% glycerol<br>pH 7.3,-20°C for 24 months(Avoid repeated<br>freeze / thaw cycles.) |
| <b>Note</b>                | For research use only                                                                                               |
| <b>Isotype</b>             | IgG                                                                                                                 |
| <b>Clonality</b>           | Polyclonal                                                                                                          |
| <b>Purity</b>              | ≥95% as determined by SDS-PAGE                                                                                      |
| <b>MW</b>                  | 60-62 kDa                                                                                                           |
| <b>Uniprot ID</b>          | <b>O76050</b>                                                                                                       |
| <b>Dilution Range</b>      | WB: 1:500-1:2000; IP: 1:200-1:2000                                                                                  |
| <b>Expiration Date</b>     | 12 months from date of receipt.                                                                                     |