

Product Datasheet

GRIN2A antibody (orb629326)

Description	Rabbit polyclonal antibody to GRIN2A
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	glutamate receptor, ionotropic, N-methyl D-aspartate 2A
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	150-160 kDa
Uniprot ID	Q12879
Dilution Range	WB: 1:500 - 1:2000; IHC: 1:50 - 1:200; IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.