

Product Datasheet

NO66 antibody (orb629340)

Description	Rabbit polyclonal antibody to NO66
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	chromosome 14 open reading frame 169
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	71 kDa
Uniprot ID	Q9H6W3
Dilution Range	WB: 1:500-1:1000
Expiration Date	12 months from date of receipt.