

Product Datasheet

NSE antibody (orb629422)

Description	Mouse monoclonal antibody to NSE
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	enolase 2(gamma, neuronal)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	47 kDa
Uniprot ID	P09104
Dilution Range	WB: 1:2000-1:20000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.