

Product Datasheet

P4HA3 antibody (orb629604)

Description	Rabbit polyclonal antibody to P4HA3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	prolyl 4-hydroxylase, alpha polypeptide III
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	30 kDa, 61 kDa
Uniprot ID	Q7Z4N8
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.