

Product Datasheet

Parkin antibody (orb629664)

Description	Rabbit polyclonal antibody to Parkin
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	Parkinson disease(autosomal recessive, juvenile) 2, parkin
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	52kDa
Uniprot ID	O60260
Dilution Range	WB: 1:500 - 1:2000; IF: 1:50 - 1:200
Expiration Date	12 months from date of receipt.