

Product Datasheet

PCBD1 antibody (orb629698)

Description	Rabbit polyclonal antibody to PCBD1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	pterin-4 alpha-carbinolamine dehydratase/dimerization cofactor of hepatocyte nuclear factor 1 alpha
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	12 kDa
Uniprot ID	P61457
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200; IF: 1:20- 1:200
Expiration Date	12 months from date of receipt.