

Product Datasheet

PDE3B antibody (orb629756)

Description	Rabbit polyclonal antibody to PDE3B
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IP
Immunogen	phosphodiesterase 3B, cGMP-inhibited
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	124 kDa
Uniprot ID	Q13370
Dilution Range	IP: 1:200-1:2000
Expiration Date	12 months from date of receipt.