

Product Datasheet

PDI antibody (orb629775)

Description	Rabbit polyclonal antibody to PDI
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, FC, IF, IHC, IP, WB
Immunogen	prolyl 4-hydroxylase, beta polypeptide
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	57kd
Uniprot ID	P07237
Dilution Range	WB: 1:500-1:1000, IP: 1:200-1:1000, IHC-P: 1:100-1:400, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.