

Product Datasheet

PDZD2 antibody (orb629798)

Description	Rabbit polyclonal antibody to PDZD2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC
Immunogen	PDZ domain containing 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	O15018
Dilution Range	IHC: 1:20-1:200; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.