



Product Datasheet

LASS2 antibody (orb6298)


Description	LASS2 antibody
Species/Host	Rabbit
Reactivity	Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC-Fr, IHC-P, WB
Immunogen	KLH conjugated synthetic peptide derived from human LASS2 (201-300/380 aa)
Target	Ceramide synthase 2
Preservatives	0.01M TBS(pH7.4) with 1% rAlbumin, 0.03% Proclin300 and 50% Glycerol.
Form/Appearance	Liquid
Concentration	1mg/ml
Storage	Shipped at 4°C. Store at -20 °C for one year. Avoid repeated freeze/thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
MW	45 kDa
Uniprot ID	Q96G23
Expiration Date	12 months from date of receipt.



Immunohistochemical staining of rat kidn...



Western blot analysis of rat kidney (Lin...



Immunohistochemical staining of rat brai...