

Product Datasheet

PILRA antibody (orb629928)

Description	Rabbit polyclonal antibody to PILRA
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA
Immunogen	paired immunoglobulin-like type 2 receptor alpha
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q9UKJ1
Expiration Date	12 months from date of receipt.