

Product Datasheet

PMS2 antibody (orb630037)

Description	Rabbit polyclonal antibody to PMS2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	PMS2 postmeiotic segregation increased 2(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	100 kDa
Uniprot ID	P54278
Dilution Range	WB: 1:500-1:2000; IHC: 1:50-1:500
Expiration Date	12 months from date of receipt.