

Product Datasheet

PPIL1 antibody (orb630124)

Description	Rabbit polyclonal antibody to PPIL1
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	peptidylprolyl isomerase(cyclophilin)-like 1
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	18 kDa
Uniprot ID	Q9Y3C6
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20- 1:200
Expiration Date	12 months from date of receipt.