

Product Datasheet

PSMB11 antibody (orb630296)

Description	Rabbit polyclonal antibody to PSMB11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	proteasome(prosome, macropain) subunit, beta type, 11
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	32-35 kDa
Uniprot ID	A5LHX3
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.