

Product Datasheet

RAB11 antibody (orb630405)

Description	Rabbit polyclonal antibody to RAB11
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	RAB11A, member RAS oncogene family
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	24 kDa
Uniprot ID	P62491
Dilution Range	WB: 1:500-1:5000, IP: 1:200-1:2000, IHC-P: 1:50-1:500, IF: 1:50-1:500
Expiration Date	12 months from date of receipt.