

Product Datasheet

RAB27A antibody (orb630423)

Description	Mouse monoclonal antibody to RAB27A
Species/Host	Mouse
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, IP, WB
Immunogen	RAB27A, member RAS oncogene family
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	25kd
Uniprot ID	P51159
Dilution Range	WB: 1:500-1:2000, IP: 1:500-1:2000, IHC-P: 1:50-1:500, IF: 1:10-1:100
Expiration Date	12 months from date of receipt.