

Product Datasheet

RAB3GAP2 antibody (orb630440)

Description	Rabbit polyclonal antibody to RAB3GAP2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	RAB3 GTPase activating protein subunit 2(non-catalytic)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	110 / 130 / 150 kDa
Uniprot ID	Q9H2M9
Dilution Range	WB: 1:500-1:2000; IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.