

Product Datasheet

RBM19 antibody (orb630565)

Description	Rabbit polyclonal antibody to RBM19
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, IHC, WB
Immunogen	RNA binding motif protein 19
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	120kd
Uniprot ID	Q9Y4C8
Dilution Range	WB: 1:500-1:5000, IHC-P: 1:20-1:200, IF: 1:10- 1:100
Expiration Date	12 months from date of receipt.