

Product Datasheet

REPS2 antibody (orb630631)

Description	Rabbit polyclonal antibody to REPS2
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, IF, WB
Immunogen	RALBP1 associated Eps domain containing 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	72 kDa, 55 kDa, 30 kDa
Uniprot ID	Q8NFH8
Dilution Range	WB: 1:500-1:2000; IF: 1:20-1:200
Expiration Date	12 months from date of receipt.