

Product Datasheet

REST antibody (orb630633)

Description	Rabbit polyclonal antibody to REST
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IP, WB
Immunogen	RE1-silencing transcription factor
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	122 kDa, 200kd
Uniprot ID	Q13127
Dilution Range	WB: 1:500-1:1000, IP: 1:500-1:2000
Expiration Date	12 months from date of receipt.