

Product Datasheet

RPE65 antibody (orb630774)

Description	Rabbit polyclonal antibody to RPE65
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	retinal pigment epithelium-specific protein 65kDa
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	65kd, 61kd
Uniprot ID	Q16518
Dilution Range	WB: 1:500-1:5000
Expiration Date	12 months from date of receipt.