

Product Datasheet

RRAS antibody (orb630861)

Description	Rabbit polyclonal antibody to RRAS
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	related RAS viral(r-ras) oncogene homolog
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	25 kDa
Uniprot ID	P10301
Dilution Range	WB: 1:500-1:3000, IP: 1:200-1:1000, IHC-P: 1:20-1:200
Expiration Date	12 months from date of receipt.