

Product Datasheet

RRN3 antibody (orb630868)

Description	Rabbit polyclonal antibody to RRN3
Species/Host	Rabbit
Reactivity	Human, Mouse
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	RRN3 RNA polymerase I transcription factor homolog(<i>S. cerevisiae</i>)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	74 kDa
Uniprot ID	Q9NYV6
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.