

Product Datasheet

S100A16 antibody (orb630921)

Description	Rabbit polyclonal antibody to S100A16
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC
Immunogen	S100 calcium binding protein A16
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
Uniprot ID	Q96FQ6
Dilution Range	IHC: 1:20-1:200
Expiration Date	12 months from date of receipt.