

Product Datasheet

SALL4 antibody (orb630945)

Description	Rabbit polyclonal antibody to SALL4
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	sal-like 4(Drosophila)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3, -20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	66-75kDa
Uniprot ID	Q9UJQ4
Dilution Range	WB: 1:500-1:2000; IP: 1:200-1:1000; IHC: 1:20- 1:200
Expiration Date	12 months from date of receipt.