

Product Datasheet

SARS antibody (orb630970)

Description This gene belongs to the class II amino-acyl tRNA family. The encoded enzyme catalyzes the transfer ...

Species/Host Rabbit

Reactivity Human, Mouse, Rat

Conjugation Unconjugated

Tested Applications ELISA, IF, IHC, WB

Immunogen seryl-tRNA synthetase

Form/Appearance liquid

Storage PBS with 0.02% sodium azide and 50% glycerol pH 7.3 , -20°C for 24 months (Avoid repeated freeze / thaw cycles.)

Note For research use only

Isotype IgG

Clonality Polyclonal


Purity ≥95% as determined by SDS-PAGE

MW 60 kDa


Uniprot ID [P49591](#)

Dilution Range WB: 1:500-1:2000, IHC-P: 1:50-1:200, IF: 1:50-1:200

Expiration Date 12 months from date of receipt.



Immunohistochemistry of paraffin-embedde...



Western blot analysis of HEK-293 cells u...