

Product Datasheet

SERPINE2 antibody (orb631081)

Description	Mouse monoclonal antibody to SERPINE2
Species/Host	Mouse
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	serpin peptidase inhibitor, clade E(nexin, plasminogen activator inhibitor type 1), member 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG1
Clonality	Monoclonal
Purity	≥95% as determined by SDS-PAGE
MW	44 kDa
Uniprot ID	P07093
Dilution Range	WB: 1:500-1:2000
Expiration Date	12 months from date of receipt.