

Product Datasheet

SIRT6 antibody (orb631202)

Description	Rabbit polyclonal antibody to SIRT6
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, WB
Immunogen	sirtuin(silent mating type information regulation 2 homolog) 6(S. cerevisiae)
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	49kd
Uniprot ID	Q8N6T7
Dilution Range	WB: 1:200-1:2000, IHC-P: 1:50-1:200
Expiration Date	12 months from date of receipt.