

Product Datasheet

SLC11A2 antibody (orb631221)

Description	Rabbit polyclonal antibody to SLC11A2
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, WB
Immunogen	solute carrier family 11(proton-coupled divalent metal ion transporters), member 2
Form/Appearance	liquid
Storage	PBS with 0.02% sodium azide and 50% glycerol pH 7.3,-20°C for 24 months(Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Purity	≥95% as determined by SDS-PAGE
MW	60-70, 75-100 kDa
Uniprot ID	P49281
Dilution Range	WB: 1:200-1:2000
Expiration Date	12 months from date of receipt.